






DISPOSABLE 3 PLY FACE MASK

-  **Food Handling**
-  **Beauty and Salon**
-  **Laboratory/General Purpose**



AVAILABLE COLORS



FEATURES & DESCRIPTION

- 3 PLY Disposable Ear Loop Face Masks provides protection.
- Secure elastic straps for easy fit.
- Pleated design for comfortable facial fit to protect from debris and particles.
- Breathable material for easy respiration.
- It is soft and comfortable, hypoallergenic without glass fibers.

Our face mask is made of 3 protective layers of high-quality non-woven fabric. The high-density filter layer can block more than 95% of particles. The inner layer provides comfort against skin and is equipt with soft elastic ear loops. We find our one size fits most users comfortably.

TECH SPECS

Product: Disposable 3 ply face mask
Material: Non-woven fabric
Color: Blue
Size: 17.5 x 9.5 cm (OSFM)
Weight: 2.8 g
Thickness: 3-ply
Standard: YY/T 0969-2013
MOQ: 2000 pcs

PACKAGING

Brand: Safe Handler
COO: PRC
Series: Everyday Series –Industrial Safety
Type: 50 pc / Box / Carton
Inner Pack Qty: 50 pc / Box
Case Pack Qty: 40 box / CTN (2000PC)
Inner Pack Size: 19x10x8.5 cm
Carton Size: 52x39x36 cm
Carton Weight: 17.64 lb
Pallet: –

PRODUCT SKU	INNER PACK UPC	CARTON GTIN
SH-D3PFM-BL-50	810025606169	
SH-D3PFM-BL-40		10810025606166

#375 Winkler Dr, Suite E, Alpharetta, GA 30004 | www.bisonlife.com | info@bisonlife.com | 1-800-201-3139

CONFIDENTIAL: This Presentation and the information contained herein is confidential and has been prepared by Bison Production Company for the purpose of setting out information relating to Bison products, business activities, patentable material or strategy. This Presentation and the information contained herein is for the exclusive use of Bison Production Company and the recipients to whom it is addressed and should not be disclosed to any third party without prior written consent from Bison Production Company.




USER INSTRUCTIONS

1. Apply the mask covering your nose, mouth and chin.
2. Apply the mask tighten over your nose.
3. Fit the mask to your face and under your chin.
4. Adjust the mask for proper wear.

⚠ WARNING

This mask does not eliminate the risk of contracting any disease or infection. Replace immediately if damaged or soiled with blood or other bodily fluids. Do not use if hypersensitive to copper or zinc. **Misuse may result in sickness or death.** Please use product as intended.



INNER BOX & OUTER CARTON



#375 Winkler Dr, Suite E, Alpharetta, GA 30004 | www.bisonlife.com | info@bisonlife.com | 1-800-201-3139

CONFIDENTIAL: This Presentation and the information contained herein is confidential and has been prepared by Bison Production Company for the purpose of setting out information relating to Bison products, business activities, patentable material or strategy. This Presentation and the information contained herein is for the exclusive use of Bison Production Company and the recipients to whom it is addressed and should not be disclosed to any third party without prior written consent from Bison Production Company.

